

Laser beam couplers 60SMS



Laser Beam Couplers 60SMS-... Lens Types and Order Codes

Lens types for Laser Beam Coupler 60SMS-...

The coupling lenses provided by Schäfter+Kirchhoff are corrected for spherical aberration and, thus, optimized for the diffraction-limited focussing of monochromatic laser radiation over a wide wavelength range. Three different kinds of optics are available, type M (laser monochromats, doublets), type A (aspheres) and type RGBV (apochromat). The focus position of type A and M optics varies with wavelength (more with aspheres than with monochromats) so that the laser beam coupler has to be refocused after any changes to the wavelength. Aspheres are designed for use with a single wavelength, while monochromats can also be used for two wavelengths over a limited range, in some cases, see Table 1, row 5.

RGBV Laser Beam Coupler 60SMS-1-4-RGBV11-47

Conventional lenses and optics exhibit chromatic aberration with a focal length dependent on wavelength (Fig. 3 [A],[B]). The coupling of combined polychromatic beams requires that chromatic aberration is corrected to the highest degree possible. For efficient coupling, all beams have to be focussed at the fiber end face onto a common point, as slight deviations lead to significant coupling losses. Any correction attempts using standard achromatic optics are really only suitable for just two different wavelengths. To circumvent these problems, Schäfter+Kirchhoff has developed their own special apochromatically corrected optics for the efficient coupling of combined polychromatic beams (RGBV, [C]) in the spectral range of 400–660 nm.

For information on the fiber-coupled Schäfter+Kirchhoff RGBV laser beam combiners, see page 51.




Figure 3:
Different types of lens for laser beam coupler 60SMS-...

Lenses offered for Laser Beam Coupler 60SMS-...

Table 1		Laser Beam Coupler 60SMS-...																				
row	curr. no.	1**	2**	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
1	Lens type	A1.5	A2	A2.7	A3.1	M3.1	M4	A4	A4.5S	M5	A6.2S	A7.5	A8	M8	A11	M11	RGBV11	M12	M12NIR	A15	M15	A18
2	Focal length f'	1.5	2	2.75	3.1	3.1	4	4	4.5	5.1	6.16	7.5	8	8.1	11	11	11	12	12	15.4	15	18.4
3	Numerical aperture NA	0.5	0.5	0.55	0.68	0.25	0.25	0.6	0.42	0.25	0.3	0.3	0.3	0.15	0.25	0.23	0.18	0.23	0.23	0.16	0.18	0.15
4	Clear aperture max. [mm]	1.5	2	3.6	5	1.7	2	5	3.7	2.5	3.7	4.5	4.9	2.5	5.5	5	4.0	5.5	5.5	5	5.5	5.5
5	Correction achromatic					x	x			x				x		x	RGBV		x		x	

Table 1.1		Input beam diameter [mm] (13.5% level)																				
6		0.04	**	**	0.22	0.25	0.32	0.36	0.40	0.50	0.60	0.64	0.64	0.88	0.88	0.88	0.88	0.96	0.96	1.23	1.20	1.47
7	Effective numerical aperture of the fiber NA _{eff} (13.5% level)	0.05	**	**	0.27	0.31	0.31	0.40	0.40	0.45	0.50	0.62	0.75	0.80	0.80	1.10	1.10	1.10	1.20	1.20	1.54	1.80
8		0.06	**	**	0.32	0.37	0.37	0.48	0.48	0.54	0.60	0.74	0.90	0.96	0.96	1.32	1.32	1.32	1.44	1.44	1.85	2.21
9		0.07	**	**	0.38	0.43	0.43	0.56	0.56	0.63	0.70	0.87	1.05	1.12	1.12	1.54	1.54	1.54	1.68	1.68	2.16	2.58
10		0.08	**	**	0.43	0.50	0.50	0.64	0.64	0.72	0.80	0.99	1.20	1.28	1.28	1.76	1.76	1.76	1.92	1.92	2.46	2.94
11		0.09	**	**	0.49	0.56	0.56	0.72	0.72	0.81	0.90	1.12	1.35	1.44	1.44	1.98	1.98	1.98	2.16	2.16	2.77	3.31
13	Multimode: nominal NA	0.22	0.66	0.88	1.19	1.36	1.36	1.76	1.76	1.98	2.20	2.73	3.30	3.52	-	4.84	4.84	-	5.28	5.28	-	-

Table 1.2		Spectral range																				
row	code no. of AR coating	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34
14	350 - 460 nm	52								52	52											
15	400 - 600 nm	01		01	01					01	01		01	01							01	01
16	600 - 1050 nm	02		02	02					02	02		02	02							02	02
17	1050 - 1550 nm	03	03	03	03					03	03		03	03							03	03
18	1300 - 1750 nm	45			45	45				45	45		45	45							45	45
19	1750 - 2150 nm	09				09				09	09		09	09							09	09
20	390 - 670 nm	33					33				33			33						33		
21	630 - 1080 nm	10								10				10					10	10		
22	980 - 1600 nm	08								08				08					08			
23	420 - 700 nm	26																			26	
24	750 - 1550 nm	37																			37	
25	400 - 670 nm	51				51												47				
26	460 - 740 nm	53												53								
27	520 - 830 nm	18														18						
28	650 - 1150 nm	07											07									
29	450 - 700 nm	04					04															
30	1750 - 3000 nm	64						64***									64***					

** for multimode fibers only, *** IR chalcogenide lens

Order Options for laser beam couplers 60SMS-...

Laser Beam Couplers 60SMS-...

60SMS-1-4-M5-33

Order Code

1 = connector type FC (standard), see dimensions 1

Optional:

AVIO = DIN-AVIO connector, see dimensions 2

ST-0-... = ST connector, see dimensions 3

SMA-... = F-SMA connector, see dimensions 4

AR coating, see Table 1.2

Lens focal length f' in mm

Lens type: A = aspheric (beam parameters Table 1)

M = laser monochromat or achromat

RGBV = laser apochromat

4 = inclined coupling axis, e.g. Fig. 2 with APC connector (8°-polish)


0 = coaxial coupling axis, e.g. Fig. 2 with PC connector (0°-polish)

Option: Add Ti for titanium construction (amagnetic): Example Order Code 60SMS-1-4-M5-33-Ti

Add SS for steel construction: Example Order Code 60SMS-1-4-M5-33-SS









Anti-Reflective Coatings for Laser Beam Couplers 60SMS-... and Fiber Collimators 60FC-...

The lenses in table 1 (page 10) and tables 1–19 (pages 35 – 49) can be ordered with the following AR-coatings (selection):



Adapters and Tools for Laser Beam Coupler 60SMS-...

Descriptions and technical drawings are available from www.SuKHamburg.com


Adapter Order Code 60A19.5-F 	Adapter Order Code 19.5 AM25 	Adapter Order Code 19.5 AM25-L 	Adapter with integrated attenuator Order Code 60A19.5-F-AT 	
Mounting Set Order Code 60A19.5-F-MS Mounting set for 60A19.5-F-... and HeNe laser 4 pcs. screws 4-40 x 3/8" (similar to DIN 912), washers and hex key 3/32 Set is supplied without adapter 	Plate Order Code 48MC-MP-19.5 For Ø 19.5 mm components, compatible with the microbench system 	Adapter with integrated shutter Order Code 60A19.5-F-S 		
Tools for assembly and adjustment <table border="0"> <tbody> <tr> <td data-bbox="119 1220 391 1377"> A Eccentric key Order Code (no. 1–16) 60EX-4 (no. 17–21) 60EX-5 as an alternative: B Eccentric key with long handle (no. 1–16) 60EX-4-L (no. 17–21) 60EX-5-L </td> <td data-bbox="399 1220 782 1377"> C Screwdriver Order Code Ø1.2 mm 9D-12 D Hex screwdriver SW Ø 1.5 mm 50HD-15 </td> </tr> </tbody> </table> 		A Eccentric key Order Code (no. 1–16) 60EX-4 (no. 17–21) 60EX-5 as an alternative: B Eccentric key with long handle (no. 1–16) 60EX-4-L (no. 17–21) 60EX-5-L	C Screwdriver Order Code Ø1.2 mm 9D-12 D Hex screwdriver SW Ø 1.5 mm 50HD-15	Iris Diaphragm   For adjustment details of laser beam coupler 60SMS-..., see page 13. Order Code 13BL1-13
A Eccentric key Order Code (no. 1–16) 60EX-4 (no. 17–21) 60EX-5 as an alternative: B Eccentric key with long handle (no. 1–16) 60EX-4-L (no. 17–21) 60EX-5-L	C Screwdriver Order Code Ø1.2 mm 9D-12 D Hex screwdriver SW Ø 1.5 mm 50HD-15			



Fundamentals: Laser Beam Couplers

Selection of lens focal length

Optimum coupling efficiency is reached when the beam diameter of the incoming beam is equal to one of the beam diameters listed in Table 1.1. Maximum coupling efficiency is achieved for an ideal Gaussian beam ($M^2 = 1$, no astigmatism) when the convergence of the focused beam equals the effective NAe^2 of the fiber [D]. Then the laser spot on the fiber end face equals the mode field diameter MFD of the single mode fiber. Except for an 8% loss from Fresnel reflection at the entrance into and



exit from the fiber, an ideal Gaussian beam is transported completely. For a specified effective fiber NAe^2 the optimum focal length of the laser beam coupler at a given beam diameter $\varnothing_{beam}$ (defined at its $1/e^2$ level) is given by

$$f' = 0.5 \cdot \varnothing_{beam} / NAe^2. \quad (1)$$

If the effective NAe^2 of the fiber is not known, then the optimum focal length f' can be calculated from the nominal numerical aperture NA by

$$f' = F_{NA} \cdot \varnothing_{beam} / NA. \quad (2)$$

The nominal fiber NA corresponds to the Gaussian angle distribution at a 1% - 5% level requiring the factor F_{NA} to correct for the different definitions of the NA .

Level	F_{NA}
1 %	0.76
3 %	0.66
5 %	0.61


Example: Beam diameter: $\varnothing_{beam} = 1.0 \text{ mm}$
Effective numerical aperture of fiber: $NAe^2 = 0.08$
focal length: $f' = 0.5 \cdot 1.0 \text{ mm} / 0.08 = 6.25 \text{ mm}$

Therefore, select e.g. lens A6.2S with $f' = 6.2 \text{ mm}$

Numerical Aperture NA

The numerical aperture NA of the optics is defined by its clear aperture. The NAe^2 of a singlemode fiber is given at its $1/e^2$ level. The NA of the coupling optics must be larger than the declared NAe^2 of the fiber. Otherwise the beam is truncated by the optics and the fiber coupling efficiency is reduced.

Selection of coupling diameter for an elliptical laser beam



In order to find the best coupling focal length in case of an elliptical beam use in eq. 1 or 2 the effective beam diameter $\varnothing_{eff}$ which is calculated from the small and the large diameters $\varnothing_{||}$ and $\varnothing_{\perp}$ of the collimated elliptical laser beam:

$$\varnothing_{eff} = \sqrt{\varnothing_{||} \cdot \varnothing_{\perp}}$$

For methods to increase the coupling efficiencies of elliptical laser diode beam sources, see page 51.

Multimode fiber


When using a multimode fiber, the coupling focal length is calculated from the beam diameter and the fiber NAe^2 :

$$f' = 0.5 \cdot \varnothing_{beam} / NAe^2.$$

A coupling focal length too long can cause mode mixing, resulting in unwanted beam characteristics, while a focal length too short will reduce the coupling efficiency.

Long term Stability of the Laser Beam Coupler

Figure 4a




The high stability of fiber-coupling using a laser beam coupler is demonstrated in temperature-stability tests using different focal lengths and wavelengths. The test setup is depicted in Fig. 4a.

The light emitted by the temperature-stabilized laser diode beam source 48TE (with integrated Faraday isolator FI) is guided to the test setup using a polarization-maintaining fiber, collimated by laser beam coupler, and then coupled back into a polarization-maintaining fiber using a second laser beam coupler, with both placed 12 mm apart. The recoupled power is monitored using a photo detector. The coupling setup is placed on a thermo-controlled plate, to vary the temperature between 15 °C and 35 °C in successive cycles with a rate of 0.5 °C per minute. The temperature of the coupling system is monitored by a temperature sensor placed on one of the two laser beam couplers.

In order to minimize any temperature impact on the measurement equipment, the laser source as well as the photo detector and the data logger are all placed on a thermo-controlled plate at a constant temperature of 25 °C.

Fig 4b shows the typical results of the relative transmitted power over 5 measurement cycles using a focal length of 4.5 mm and a wavelength of 405 nm. The power is normalized with respect to the mean power acquired over all measurement cycles. The power deviation from the mean power is $\pm 1.5 \%$.

The repetitive pattern in the relative power caused by the temperature



cycling is demonstrated more clearly in Fig. 4c, in which the relative power (normalized to the maximum) is plotted against the temperature of the laser beam couplers. In this case the maximum coupling efficiency is reached a little above 25 °C and it decreases faster towards lower temperatures than higher temperatures, with the smallest slope near the requested operating point (25 °C).

The respective power curves for each measurement cycle are almost coincident and the power variation at points with equal temperatures is $< 1 \%$, which demonstrates the reproducibility of the pointing stability during temperature cycling and the long-term stability of the fiber-coupling.

The maximum deviation with respect to the maximum power here is 3 %.




Figure 4:

Test setup for measuring the stability of two laser beam couplers ($f = 4.5 \text{ mm}$, $\lambda = 405 \text{ nm}$) during successive temperature cycling between 15 °C and 35 °C (a). The relative power (normalized with respect to the max power) (b) shows a repetitive pattern following the temperature (below) and has a maximum deviation of $\pm 1.5 \%$.

The relative power curves (c) (normalized with respect to the maximum power) are almost coincident and confirm the high reproducibility of the pointing stability during temperature cycling. The maximum deviation is only 3 %.


Laser Beam Coupler 60SMS-... Assembly and adjustment

For efficient beam coupling, the choice of the coupling optics is made according to the diameter of the incoming beam and the numerical aperture of the coupled fiber. The precision adjustment mechanism is then used for the precise placement of the laser spot on the transmitted mode field of the fiber end face and – with polarization-maintaining fibers – for orientation of the polarization-maintaining axis of the fiber so that it is fully parallel with the polarization plane of the incoming radiation.


Four steps for perfect alignment:

1. Center the laser beam coupler with the laser beam propagation axis **[1]** by using the adapter 60A19.5-F (or similar).
2. Move the mode field of the fiber laterally within the laser focus area using the tilt adjustment **[2]**.
3. Alter the preset factory-setting of the focal plane **[3]** (only needed for a wavelength different from specification or when beam splitting is experienced).
4. Rotate the laser beam coupler to align the polarization-maintaining axes (only PM-fibers) **[4]**.

Centering of laser beam coupler with the propagation axis **[1]**



Coupling efficiencies with f' 5–15 m lenses for uncentered parallel beams




(using the deliberately oversized mounting holes) while measuring the laser power.

The centering precision required is of the order of tenths of a millimeter and, with a longer focal length lens, manual adjustment is perfectly adequate and quickly performed **[A]**.

Lateral adjustment of the mode field and laser spot **[2]**

Lateral displacement of the laser beam focus from the mode field of the fiber arises because of:

- production tolerances in the centering of the coupling lens and/or the centering of the fiber core in the fiber ferrule. With a mode field or spot diameter of 2–5 μm , the required precision is in the sub-micron range **[D]**.



Coupling efficiencies for laser spots laterally displaced from the fiber mode field

- inclined beam propagation **[E]**.

Using a 5 mm focal length lens, a beam inclined by 1 mrad results in a lateral offset of 5 μm – completely


missing the mode field. By using the tilt mechanism of the laser beam coupler, the mode field of the fiber is then placed very precisely onto the laser focus area.

Cyclical tightening and release **[F]** of the three radially arranged adjustment screws

maximize the coupled beam power. The final position is locked by the adjacent grub screws, which have carbide bearings.

Focussing **[3]**

The positioning accuracy in the coaxial direction is less critical than for the lateral directions. From the depth of field of the laser focus (Rayleigh range), a decrease in coupling efficiency of a few percent occurs even with a displacement of only a few microns **[G]**.



Coupling efficiency with unfocused lens




eccentric key 60EX-4 to relocate the focal plane, before tightening the locking screws again.

Schäfter + Kirchhoff laser beam couplers are supplied adjusted for the coupled wavelength and refocussing is not necessary for a highly collimated laser beam.

Refocussing can be done, however, by releasing the two lens-locking screws (accessible via small holes) using screwdriver 9D-12 **[H]**. Use the

1 Cross-section of a laser beam coupler and its functional and adjustable elements

- 1 Laser beam coupler 60SMS-1-4-...
- 2 Adapter 60A19.5-F
- 3 Laser beam source (HeNe laser)
- 4 Singlemode fiber cable with FC-APC connector
- 5 Grub screws for fixing the fiber ferrule position
- 6 Coupling lens
- 7 Lens focussing mechanism
- 8 Locking screws for the indirect fixing of the coupling lens
- 9 Adjustment screws (3 x tilt adjustment)
- 10 Locking screws (3 x tilt adjustment)
- 11 Carbide bearings (3 x)





Orientation of the polarization axis **[4]**

Polarization-maintaining singlemode fibers guide radiation in the two orthogonal planes of oscillation (fast and slow), independently. When linearly polarized radiation is not orientated exactly with one of the fiber axes, the radiation is split across the fast and slow propagation speeds and recombines at the fiber end as elliptically polarized radiation; the linear state of polarization is lost **[J]**.

The polarization state of a fiber can be influenced by bending, vibration and temperature.

A linear and stable state of polarization is obtained by rotating the laser beam coupler in the adapter 60A19.5-F **[K]** (after releasing the locking screws), and confirmed by the lack of effect of gentle physical displacement of the fiber. This new position can now be locked by retightening the screws.



Extinction ratio for polarization axis twisted in respect of the slow and fast axes of the fiber

for analysis of the state of polarization at the fiber end. The analyzer evaluates the extinction ratio and immediately displays the results of any adjustments made to the polarization.

The SK010PA-... polarization analyzer was specifically developed



Polarization Analyzer SK010PA-...

For measuring and testing laser beams with polarization-maintaining singlemode fibers, see page 71